D

Brilliance®

Gebrauchsinformation

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Brilliance®

Instructions for use

Addition-curing bite registration material based on polyvinyl siloxane. Very high final hardness.

Indications and properties

- Syringeable bite registration material, which is also indicated for use in centric occlusion.
- Material for coating bite registration forks.
- A locating material, such as for Gothic arch tracings.

Thanks to its high final hardness, this material is easily milled. Thanks to its thixotropic properties, Brilliance is easy to apply. Its firm consistency prevents it flowing into the interdental spaces and enables it to be withdrawn easily from the mouth. We recommend drying the site prior to application. Brilliance bite registrations exhibit permanent dimensional stability.

Instructions for use with the dispensing gun

- 1) Insert the cartridge into the dispensing gun. Ensure that the notches on the cartridge base are pointing downwards. The clip will not close if the cartridge is not properly inserted into the delivery system.
- 2) Remove the cartridge closure.
- 3) Do not mount the mixing tip. Extrude impression material until it exits both chambers at the same rate.
- 4) Mount the mixing tip on the cartridge and extrude the required amount of material.
- 5) Leave the filled mixing tip on the cartridge after use or apply the cartridge closure.
- 6) Before using the cartridge again, remove and discard the used mixing tip/the cartridge closure. Check the cartridge openings to ensure that no plugs of hardened material remain. Should this be the case, unblock the opening the openings with a pointed instrument and check that the material is extruded from both at the same rate.
- 7) Mount a fresh mixing tip and proceed as usual.

The dispensing gun can be sterilized in the autoclave.

Technical Data

Total working time:	a minimum of 30 sec.
Intraoral setting time:	a minimum of 1 min.
Setting time:	2 min.
Shore-A hardness:	90
Linear dimensional change:	-0,20%

Storage conditions:

In a dry place at 18-25°C (64-77°F).

Disinfection:

Bite registrations can be disinfected.

Expiry date:

See cartridge

To ensure excellent bite registrations, the total working time and intraoral setting time must be adhered to exactly. Should the temperature of the material deviate noticeably from 23°C (74°F), these times will be affected.

Warning

Solutions applied prior to taking bite registrations (e.g. mouthwashes), may impede the setting reaction of the registration material and must be cleaned off by spraying thoroughly with water and drying.

Do not leave any remains of the registration material in the oral cavity.

Avoid splattering this material onto clothing as it cannot be removed by dry cleaning.

The possibility of allergic reactions to the product in sensitive persons cannot be fully excluded.

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Willmann & Pein GmbH