HR

## Securafil®

### Upustvo za upotrebu

Securafil je stakloionomeri cementni materijal za punjenje dizajniran za upotrebu u klas III i V šupljina, kao što su erozije u obliku slova V, karjes u korijenu, karjes na ivicama krunika, karjes na aproksimalnim tačkama kontakta (pripremanje tunela) i u zubnom mesu. Securafil se može takođe koristiti kod zaščitovanja zubnih pukotina. Securafil se veže sa dentinom i caklinom pomoću reakcije karboksilne grupe polikarbonske kiseline sa kalcijem hidroksipatita i kolagenom. Ova hemijska reakcija čuva mnogo više zubnog materijala tako da bilo kakva podređivanja zuba pri dnu nisu neophodna. Stakloionomeri punjenje sadrži fluorid.

Securafil je dostupan u sljedećim nijansama boje: A2, A3

### Priprema kanala

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## Securafil®

### Instructions for use

Securafil je a glass ionomer cement filling material designed for use in class III and V cavities, such as V-shaped erosions, root caries, caries in crown margins, caries in the approximal contact points (tunnel preparation), and in delicious teeth. It may also be used as fissure sealant. Securafil bonds to dentine and enamel by reaction of the carboxylic acid with hydroxyapatite calcium and with collagen. This chemical reaction preserves sound material. No undercuts are necessary. The glass ionomer filler contains fluoride.

Securafil is available in following colourshades: A2, A3.

**Side effects**

Undesirable side-effects to this medical product are not to expect, if working corresponding to the instructions of use. In case of side-effects - even when doubtful - please inform us with exact description circumstances and symptoms.

### Instruction references

Take care of all working- and caution advices. In case of unexpected side effects or other reactions (even if work was made according to the instructions of use) please inform us, if possible, with exact description of the accompanying circumstances and symptoms.

### Cavity preparation

Sound dentine should be preserved as far as possible since undercuts are not necessary due to Securafil's good adhesion. Cavity design preferably should allow the edge thickness of the cement layer to be more than 0,5 mm.

### Pulp protection

Cavities are rinsed carefully with water and dried with oil free air. Deep cavities should be protected with a calciumhydroxide liner.

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### Application of cleanser

Apply tooth cleanser on freshly prepared dentine surfaces or use it to remove plaque. Leave the tooth cleanser on the surfaces for 10 seconds, rinse with water and dry using oil free air. Avoid Dehydration!

### Mixing of Securafil

Shake powder before dispensing. Use two level measuring scoops of powder for two full drops of liquid (without air-bubbles, hold bottle upright, press bottle slightly). Mix powder into liquid beginning with half of the powder and continue taking half of the remaining powder and so on. To minimize evaporation of the liquid use only a small area on the mixing pad.

### Application of Securafil

The ready mixed material is applied into the cavity. The material is condensed slightly. The surface is shaped using a matrix. If used without matrix, it should be covered with an ionomer varnish. Securafil is malleable for 2 min (incl. Mixing). Once the hardening process starts, the cement should not be disturbed. Setting time is 3 min.

### Finishing

After the setting of the matrix can be removed. The now unprotected area should be covered with an ionomer varnish. It is recommended to finish the material after additional 6 min. Excess of material can be removed with white stones or flexible discs which have been lubricated with petroleum jelly.

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